MARGARET GOMEZ
NOTARY PUBLIC
WAKE COUNTY, N.C.
My Commission Expires 9/3/27

RADIESSE LIDOCAINE

INJECTABLE IMPLANT

INSTRUCTIONS FOR USE

NAME

RADIESSE Lidocaine injectable implant,

Containing:

- injectable implant in a 1.5 ml syringe – 1 pc.
- needle for injections single-use sterile K-Pack II Needle 27G × ¾" TW (0,4 × 20 mm) (by TERUMO EUROPE N.V., Belgium, RU № RZN 2018/7086) – 2 pcs.
- instruction for use – 1 pc.

INTENDED USE

RADIESSE Lidocaine injectable implant is intended for plastic and reconstructive surgery, including deep dermal and sub-dermal soft tissue augmentation of the facial area.

INTENDED USER

RADIESSE Lidocaine injectable implant should be used by healthcare practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.

RADIESSE Lidocaine injectable implant is designed to be used in medical and preventive treatment institutions.

DESCRIPTION

RADIESSE Lidocaine injectable implant is a sterile, latex-free, semi-solid, cohesive implant. The principal component is synthetic calcium hydroxylapatite suspended in a gel carrier that consists of sterile water for injection (USP), glycerin (USP), sodium carboxymethylcellulose (USP), monobasic sodium phosphate dihydrate (USP, Ph. Eur., BP), dibasic sodium phosphate anhydrous (USP, Ph. Eur., BP) and 0.3% lidocaine hydrochloride (USP). The gel is dissipated *in vivo* and replaced with soft tissue growth, while the calcium hydroxylapatite remains at the site of injection. The period of biodegradation of RADIESSE Lidocaine implant at the site of injection is up to 30 months after injection (implantation). The lidocaine provides short-term local anesthetic effect. The result is long-term yet non-permanent restoration and augmentation.

RADIESSE Lidocaine injectable implant is manufactured in a 1.5 ml syringe, has a calcium hydroxylapatite particle size range of 25-45 microns and is injected with a K-Pack II Needle 27G ¾" TW (0,4 mm x 20 mm) needle with a standard Luer/Luer lock fitting.

PRESENTATION

RADIESSE Lidocaine injectable implant is presented in a pre-filled syringe for single use, packaged in a foil pouch for terminal steam sterilization. Each box contains one instruction for use, one 1.5 ml syringe in a foil pouch and two sterile needles 27G × ¾" TW (0,4 × 20 mm) for single use with CE-mark.

Two traceability labels are incorporated into the pouch label.

IMPLANT COMPOSITION

Calcium hydroxylapatite ($\text{Ca}_{10}(\text{PO}_4)_6(\text{OH})_2$) – 937,10 mg/ml
Glycerin ($\text{C}_3\text{H}_8\text{O}_3$) – 105,80 mg/ml
Sodium carboxymethylcellulose ($\text{C}_8\text{H}_{16}\text{NaO}_8$) – 17,60 mg/ml
Lidocaine hydrochloride ($\text{C}_{14}\text{H}_{23}\text{N}_2\text{OCl}$) – 3,70 mg/ml
Sterile water for injection ($\text{C}_{14}\text{H}_{23}\text{N}_2\text{OCl}$) – 590,80 mg/ml
Dibasic sodium phosphate anhydrous (Na_2HPO_4) – 4,50 mg/ml
Monobasic sodium phosphate dihydrate ($\text{NaH}_2\text{PO}_4 \times 2\text{H}_2\text{O}$) – 0,55 mg/ml

INFORMATION ABOUT THE PRESENCE IN THE COMPOSITION OF THE MEDICAL DEVICE OF A MEDICINAL PRODUCT FOR MEDICAL USE, MATERIALS OF ANIMAL AND OR HUMAN ORIGIN

Lidocaine hydrochloride produced by Albemarle Corporation is a medicinal product that is part of the Radiesse injection implant with lidocaine.

RADIESSE Lidocaine injectable implant does not contain any animal or human derived components.

INDICATIONS

RADIESSE Lidocaine injectable implant is indicated for

- the treatment of nasolabial folds
- the augmentation of cheeks
- the treatment of marionette lines
- the treatment of jawline

TARGET TREATED POPULATION

Adult patients regardless of gender of all ethnicities and all Fitzpatrick skin types with respect to indications and contraindications stated in the instructions for use.

CONTRAINDICATIONS

- RADIESSE Lidocaine injectable implant is contraindicated in the presence of acute and/or chronic inflammation or infection when these involve the area to be treated.
- Contraindicated in patients with severe allergies manifested by a history of anaphylaxis, or history or presence of multiple severe allergies.
- Not to be used in patients with known hypersensitivity to any of the components.
- RADIESSE Lidocaine injectable implant is not intended to be used in patients with known hypersensitivity to lidocaine or anesthetics of the amide type.

- Radiesse Lidocaine injectable implant is contraindicated in patients prone to developing inflammatory skin conditions or those patients with a tendency for developing hypertrophic scars and keloids.
- Do not implant in the epidermis or use as a skin replacement. Implantation into the epidermis or superficial dermis could lead to complications such as fistula formation, infections, extrusions, nodule formation and induration.
- Radiesse Lidocaine injectable implant is not intended to be used for the correction of glabellar folds and nose area. A higher incidence of localized necrosis has been associated with glabellar and nose injections. Complications indicate that forceful injection into superficial dermal vessels of the glabellar or nose area could cause retrograde movement into the retinal arteries resulting in vascular occlusion.
- Radiesse Lidocaine injectable implant is contraindicated in the presence of foreign bodies such as liquid silicone or other particulate materials, such as other filler containing particles are meant, e. g. HA technology products containing particles, or other CaHA particle filler.
- Radiesse Lidocaine injectable implant should not be used in areas where there is inadequate coverage of healthy, well vascularized tissue.
- Radiesse Lidocaine injectable implant should not be used in patients with systemic disorders, for example systemic lupus erythematosus, rheumatoid arthritis and scleroderma, which cause poor wound healing or will lead to tissue deterioration over the implant.
- Radiesse Lidocaine injectable implant is contraindicated in patients with bleeding disorders.

WARNINGS

- Introduction of Radiesse Lidocaine injectable implant into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction. Take extra care when injecting soft tissue fillers, for example inject the Radiesse Lidocaine injectable implant slowly and apply the least amount of pressure necessary. Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate healthcare practitioner specialist should an intravascular injection occur.
- Do not overcorrect (overfill) the injection site as the soft tissue volume is expected to increase within several weeks as the treatment effect of Radiesse Lidocaine injectable implant occurs.
- Implant should not be implanted in patients while the patient is on an aspirin regimen or while taking other medications that could inhibit the healing process.
- Implant should not be implanted in infected or potentially infected tissue or in open cavities because infection or extrusion may occur. A significant infection may result in damage or loss to the skin overlying the implant. Hematomas or seromas may require surgical drainage.
- Use of Radiesse Lidocaine injectable implant in patients with active skin inflammation or infection in or near the treatment area should be deferred until the inflammatory or infectious process has been controlled.
- In the event of a hypersensitivity or allergic reaction, a significant inflammation or infection may occur requiring the removal of the implant.

- Injectable implants have been associated with hardening of the tissues at an injection site, migration of particles from an injection site to other parts of the body and/or allergic or autoimmune reactions.
- As with any implant material, possible adverse reactions which may occur include, but are not limited to, the following: inflammation, infection, fistula formation, extrusion, hematoma, seroma, induration formation, inadequate healing, skin discoloration and inadequate or excessive augmentation.
- Safety and performance during pregnancy, in breastfeeding females, or in patients under 18 years old has not been established.

PRECAUTIONS

- In order to minimize the risks of potential complications, RADIESSE Lidocaine injectable implant should only be used by healthcare practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.
- In order to minimize the risks of potential complications, healthcare practitioners should fully familiarize themselves with the product, the product educational materials and the entire package insert.
- The safety and effectiveness of a cannula used for deep injection (subdermal and/or supraperiosteal) for soft tissue augmentation to improve the contour of the jawline has only been clinically evaluated in the following blunt-tip cannula (Sterimedix Silkann) that was 27G and 40mm in length and a 25G Pre- Hole Needle (clinical US study with RADIESSE Lidocaine Injectable Implant).
- Healthcare practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.
- As with all transcutaneous procedures, RADIESSE Lidocaine injectable implant injection carries a risk of infection. Infection may necessitate attempted surgical removal of RADIESSE Lidocaine injectable implant. Standard precautions associated with injectable materials should be followed.
- Patients who are using medications that can prolong bleeding, such as aspirin or warfarin, may, as with any injection, experience increased bruising or bleeding at the injection site.
- The calcium hydroxylapatite (CaHA) particles in RADIESSE Lidocaine injectable implant are clearly visible on CT Scans or mammograms and may be visible in standard, plain radiography. Patients need to be informed of the radio-opaque nature of RADIESSE Lidocaine injectable implant, so that they can inform their primary care health professionals and/or radiologists. In a radiographic study of 58 patients, there was no indication of RADIESSE injectable implant (without lidocaine) potentially masking abnormal tissues or was interpreted as tumors in CT Scans.
- RADIESSE Lidocaine injectable implant requires soft tissue for easy percutaneous injection. Scar tissue and significantly compromised tissue may not accept the implant appropriately.
- Infection requiring treatment may occur at the injection site. If such infection cannot be corrected, it may become necessary to remove the implant surgically.
- Injection related reactions, including bruising, erythema, swelling, pain, itching, discoloration or tenderness, may occur at the site of the injection. These usually resolve spontaneously within one to two days after the injection.

- Nodule(s) may form requiring treatment or removal.
- Irregularity of the implant may occur which may require a surgical procedure to correct.
- Do not over-inject the area to be treated. In extreme cases site rupture could occur. Radiesse Lidocaine injectable implant can be easily added in subsequent injections, and can be removed only surgically.
- The Radiesse Lidocaine injectable implant injection procedure, like similar injection procedures, has small but inherent risks of infection and/or bleeding. The usual precautions associated with percutaneous injection procedures should be followed to prevent infection.
- **Do not re-sterilize.** Radiesse Lidocaine injectable implant is supplied sterile and non-pyrogenic in a sealed foil pouch and is intended for single patient, single treatment use only.

The foil pouch should be carefully examined to verify that neither the pouch nor the syringe has been damaged during shipment. Do not use if the foil pouch is compromised or the syringe has been damaged. Do not use if the syringe end cap or syringe plunger is not in place. There is a small amount of moisture normally present inside the foil pouch for sterilization purposes; this is not an indication of a defective product.

- To help avoid needle breakage, do not attempt to straighten a bent needle. Discard it and complete the procedure with a replacement needle.
- Do not recap used needles. Recapping by hand is a hazardous practice and should be avoided.
- The safety of Radiesse Lidocaine injectable implant with concomitant dermal therapies such as epilation, UV irradiation, radiofrequency, ablative or non-ablative laser, mechanical or chemical peeling procedures has not been evaluated in controlled clinical trials.
- No studies of interactions of Radiesse Lidocaine injectable implant with drugs or other substances or implants have been conducted.
- Limited clinical data is available on the combination of Radiesse Lidocaine Injectable Implant with BELOTERO products and/or botulinum toxin. As a precaution, products should be injected in different facial areas. HCPs should be experienced, and persons appropriately selected as not only benefits, but also adverse events can be cumulative, and causality of adverse events could become difficult to determine. Instructions for use, depth of injection and appropriate recommendation of each product should be followed.
- Universal precautions must be observed when there is a potential for contact with patient body fluids. The injection session must be conducted with aseptic technique.
- If laser treatment, chemical peeling, or any other procedure based on active dermal response is considered after treatment with Radiesse Lidocaine injectable implant, there is a possible risk of eliciting an inflammatory reaction at the implant site. This also applies if Radiesse Lidocaine injectable implant is administered before the skin has healed completely after such a procedure. When combining dermal fillers with other aesthetic techniques such as peeling, dermabrasion or any type of laser treatment, procedures should be spaced by 2 weeks to allow complete healing of the last treatment.
- Injection of Radiesse Lidocaine injectable implant into patients with a history of previous herpetic eruption may be associated with reactivation of the herpes virus.
- Safety of Radiesse Lidocaine injectable implant beyond 1.5 years in the face has not been investigated in clinical trials. Care should be taken to determine the risk versus the benefit for patients with congenital methemoglobinemia with glucose-6-phosphate dehydrogenase deficiencies, and with patients who are receiving concomitant treatment with methemoglobin-inducing agents.

ADVERSE EVENTS

The patients must be informed that as with any implant material possible adverse reactions may occur.

Adverse events seen in a clinical trial with RADIESSE Lidocaine injectable implant were generally expected, mild in nature, and short in duration.

The following adverse events were reported during clinical trials performed with the RADIESSE injectable implant (without lidocaine): ecchymosis, edema, erythema, nodule, pain, pruritus, soreness, tenderness, numbness, contour irregularity, lumps, irritation, rash, needle jamming, discoloration, hardness, headache, scab, tightness, blood shot eyes, black eye, abrasion, spot, nerve sensitivity, dry, burning sensation, warm, stretched, pimple, flushed, feverish, ear running, backed-up salivary gland, firmness, hearing loss, and puffiness.

The following adverse events have been identified during post-approval use of RADIESSE injectable implant. (without lidocaine). Because they are reported voluntarily from a population (including from literature) of uncertain size, it is not always possible to reliably estimate their frequency or establish a causal relationship to RADIESSE injectable implant. (without lidocaine) These events have been chosen for inclusion due to the combination of their seriousness, frequency of reporting, or potential causal connection to RADIESSE injectable implant (without lidocaine):

infection (incl. biofilm formation), cellulitis, impetigo, loss of effect, product displacement/migration, allergic reaction, anaphylaxis, hives, rash, pruritus, urticaria, angioedema, inflammation, necrosis, granuloma, nodules, induration, erythema, skin discoloration (including hypo- and hyperpigmentation), dissatisfaction, pustule, skin pallor, hair loss, paresthesia, ptosis, pain (incl. mastication pain), headache, swelling, asymmetry, abscess, herpetic infection including herpes simplex and herpes zoster, hematoma, blanching, blistering, dizziness, festoons, flu-like symptoms, Guillain-Barre syndrome, tachypnea, ischemic reaction, lymphoid hyperplasia, lymphatic obstruction, nausea, pericarditis, scarring, sensitivity to cold, vascular occlusion/obstruction, vascular compromise, ocular ischemia, diplopia, visual impairment/blindness, facial muscle paralysis, Bell's palsy, syncope, chewing problems, erosion, injection site cyst, injection site warmth, aggravation of preexisting conditions, superficial vein prominence, vasculitis, nerve injury, nerve compression, xanthelasma, anxiety, laceration, purpura, calcinosis, capillary fragility, contact dermatitis, dysphagia, dyspnoea, ear pain, eye pain, eyelid function disorder, hypoesthesia, nerve damage, vesicles, lip disorder, mouth ulceration, muscle twitching, nasal congestion, nasopharyngitis, neoplasm, rhinorrhea, exfoliation, fissures, wounds, sneezing, toothache, vomiting.

Patients with specific ethnic characteristics, e.g. Asian population, should be informed of a higher risk of tissue reactions, e.g. inflammatory reactions, pigmentary disorders, post-inflammatory hyperpigmentation (PIH), scarring, and keloid formation upon cutaneous injury.

The following interventions have been reported: antibiotics, anti-inflammatories, corticosteroids, antihistamines, analgesics, massage, warm compress, excision, drainage, and surgery. This information does not constitute and is not intended to be medical advice, a recommendation on how to treat an adverse event or an exhaustive list of possible interventions. Physicians should evaluate each case on an individual basis, and independently determine, based on their professional experience, what treatment(s) are appropriate, if any, for their patients.

CIRCUMSTANCES IN WHICH THE PATIENT SHOULD CONSULT WITH A HEALTH CARE PROFESSIONAL

If you experience any adverse events listed in the "Adverse events" section, you should consult a health care professional.

INFORMATION ABOUT THE METHODS BY WHICH IT IS POSSIBLE TO DETERMINE THE LOCALIZATION OF THE IMPLANT

Visualization of a medical device *in vivo* is possible by means of MRI and / or ultrasound diagnostics.

INDIVIDUALIZATION OF TREATMENT

Before treatment, the patient's suitability for the treatment should be assessed. The outcome of treatment will vary among patients. In some instances additional treatments may be necessary depending on the size of the defect and the needs of the patient. Additional injections may be performed, but only after sufficient time has passed to evaluate the patient. The patient should not be re-injected sooner than seven days after the previous treatment.

DIRECTIONS FOR USE

GENERAL

For the percutaneous injection procedure the sterile needles 27G × ¾" TW (0,4 × 20 mm) with Luer/Luer lock fitting that are part of the product presentation should be used together with a 1,5 ml syringe.

1. Prepare patient for percutaneous injection using standard methods. The treatment injection site should be marked by a surgical marker and prepared with a suitable antiseptic.
2. Prepare the syringes and the injection needle(s) before the percutaneous injection. A new injection needle may be used for each syringe, or the same injection needle may be connected to each new syringe for same patient treatment.
3. Remove foil pouch from the carton. The pouch can be opened and the syringe dropped onto the sterile field when required. *There is a small amount of moisture normally present inside the foil pouch for sterilization purposes; this is not an indication of a defective product.*
4. Peel or twist apart the needle packaging to expose the hub.
5. Remove the Luer lock syringe cap from the distal end of the syringe prior to attaching the needle. The syringe can then be twisted onto the Luer lock fitting of the needle taking care not to contaminate the needle. Discard needle package. **The needle must be tightened securely to the syringe and primed with RADIESSE Lidocaine injectable implant.** If excess implant is on the surface of the Luer lock fittings, it will need to be wiped clean with sterile gauze. Slowly push the syringe plunger until the implant material extrudes from the end of the needle. If leakage is noted at the Luer lock fitting, it may be necessary to tighten the needle or to remove the needle and clean the surfaces of the Luer lock fitting or, in extreme cases, replace both the syringe and the needle.
6. Locate the initial site for the implant. Scar tissue and cartilage may be difficult or impossible to inject. Examples for scar tissue in the face are: acne scars of the face, scars from burns, dog bites, accidents or surgeries. An example for cartilage tissue in the face is the temporomandibular joint. Avoid passing through these tissue types when advancing the injection needle.

NOTE: Do not inject into a blood vessel.

7. The depth of the injection and the amount injected will vary depending on the site and extent of the restoration or augmentation. RADIESSE Lidocaine injectable implant should be injected sufficiently deep to prevent nodular formation at the surface of the skin or ischemia of the overlying tissue.
8. **DO NOT OVERCORRECT THE INJECTION SITE.** Use a 1:1 correction factor. Mold or massage the injected implant periodically during the injection process to maintain a smooth contour of the implant. A maximum of 10ml RADIESSE Lidocaine injectable implant is recommended as yearly dose. The dose may be adapted or increased according to the person's indication(s), tissue, age, injection depth and technique for implantation.
9. If significant resistance is encountered when pushing the plunger, the injection needle may be moved slightly to allow easier placement of the material. If significant resistance is still encountered, it may be necessary to pull the needle entirely out of the injection site and try again in a new position. If significant resistance continues to persist, it may be necessary to try a different injection needle. If this is not successful, replace the syringe and injection needle.
10. Advance the needle bevel down at approximately a 30° angle to the skin into the sub-dermis to the starting location. [Refer to the section "TREATMENT OF FACIAL INDICATIONS"]. Carefully push the plunger of the syringe to start the injection and slowly inject the implant material while withdrawing the needle, placing a line of material in the desired location. Continue placing additional lines of material until the desired level of augmentation is achieved. The thread of implant material should be completely surrounded by soft tissue without leaving globular deposits. The injected area may be massaged as needed to achieve even distribution of the implant.

TREATMENT OF FACIAL INDICATIONS

1. Insert needle with bevel down at approximately a 30° angle to the skin. The needle should slide into the deep dermis to the point where you wish to begin the injection. This should be easily palpable with the non-dominant hand.
2. Apply slow continuous even pressure to the syringe plunger to inject the implant as you withdraw the needle, leaving behind a single thin thread or strand of implant material. The thread of implant material should be completely surrounded by soft tissue without leaving globular deposits.
3. Individual threads of implant material should be placed parallel and adjacent to each other, and layered when deeper folds are corrected. As an option, the threads can be cross layered in a deeper plane for structural support.
4. After injection, use the index finger and thumb to smooth the areas and better distribute the implant in case of any slight nodular deposition of material.

PATIENT COUNSELING INFORMATION

The patient should be instructed in appropriate post-procedural care, which may include the following, to promote normal healing and avoid complications.

- Apply cool compresses to areas of injection for approximately 24 hours.
- Avoid the sun, tanning (ultraviolet) lights, sauna and intense facial treatments post procedurally.
- Massage the area gently if palpable nodules become present.
- Promote facial rest for one week by encouraging patients to limit talking, smiling and laughing.

- Inform the patient that postoperative swelling and numbness are common. Swelling will usually resolve within 7 to 10 days, but may persist for several weeks. Numbness should resolve within 4 to 6 weeks.

Do not use if packaging and/or syringe are damaged or if the syringe end cap or syringe plunger is not intact.

The contents of the syringe are intended for single patient, single treatment use only and cannot be re-sterilized. Re-use may compromise the functional properties of the device and/or lead to device failure. Re-use may also create a risk of contamination of the device and/or cause patient infection or cross-infection including but not limited to transmission of infectious disease(s) and blood transfer between patients. All of which, in turn, may lead to patient injury, illness or death.

STERILITY

Pre-filled syringe for single use, packaged in a foil pouch, is to be terminally sterilized by steam.

Foil pouch allows it to withstand the sterilization process and attain the required conditions for sterilization within the packaging system.

Foil pouch integrity may be used to demonstrate maintenance of the implant sterility.

Foil pouch maintains sterility of the medical device until the expiry date.

The needles are sterilized using ethylene oxide.

TECHNICAL CHARACTERISTICS

RADIESSE Lidocaine Injectable Implant specification:

| Test/Property | Values |
|---|---|
| Pouch Seal Strength | ≥5,0 lbf (≥22,24 N) |
| Pouches/Seals | 1/4" minimum seal width, no holes, voids or wrinkles that affect seal integrity |
| Media Fill | 1.5cc volume; no voids or absence of color uniformity |
| Labeling | Labels correct and properly applied, labels clear and legible |
| pH value | 6.8 – 7.6 |
| Media Percent solids | 54.3% - 58.3% (w/w) |
| Lidocaine HCl concentration | 0.20%-0.40% |
| 2,6-Dimethylaniline (2,6-DMA) concentration | ≤1.65 ppm |
| Sterility | Sterile |
| Bacterial endotoxins | <0.5 EU/ml |
| Dynamic viscosity (η) (0.4 Hz) | 219.679 Pa·s – 425.000 Pa·s (219679 cP – 425000 cP) |
| Tan Delta (0.6 Hz) | 0.46 – 0.76 |
| Final Packaging/Release | Compliant |

Molecular mass of calcium hydroxylapatite is 502 Da.

Volume of implant in each syringe is from 1,55 ml to 1,65 ml. Deliverable volume is not less than 1,5 ml.

RADIESSE Lidocaine injectable implant is presented in a pre-filled syringe for single use, packaged in a foil pouch for steam sterilization.

Surfaces of the pouch are joined together by thermal fusion. The seal provides a barrier to microorganisms during sterilization, transportation and storage under specified conditions. The seal width is not less than 6,35 mm.

The delivery package includes 2 sterile single-use thin-walled needles for injections K-Pack II Needle 27G x 3/4" TW (0,4 x 20 mm).

Since April 26, 2018, K-Pack II Needle needles have been put into circulation in the territory of the Russian Federation (certificate of registration No. RZN 2018/7086).

Dimensions of the needles supplied with RADIESSE Lidocaine injectable implant – needles for injections single-use sterile K-Pack II Needle

| Product code/ size, mm | External diameter, mm | Internal diameter, mm | Bevel angle (α), ° | Effective length of the needle cannula, mm |
|---------------------------|-----------------------|-----------------------|-----------------------------|--|
| KN2719-RBT / 0.40 x 20 mm | 0.400 – 0.420 | Not less than 0,241 | 12° (long bevel) | 20 +1/-2 |

THE LIST AND DESCRIPTION OF MATERIALS, COMING INTO DIRECT OR INDIRECT CONTACT WITH THE PATIENT'S BODY

Materials of a syringe:

| Component | Material | Direct contact with implantat |
|----------------|--|-------------------------------|
| Syringe barrel | Cyclic Olefin Copolymer, transparent No colorant | Yes |
| Syringe Cap | Cyclic Olefin Copolymer, semitransparent No colorant | Yes |
| Push Rod | Polycarbonate (10% Filled), white No colorant | No |
| Plunger | Silicone rubber Colorants: 0,25% Blue silicon paste 0,1% White silicon paste | Yes |

Syringe graduations in printed with black ink ITW Type P RNT. The degree of accuracy of syringe graduations is ± 0.025 ml.

Needle materials:

| Needle component | Material |
|------------------|--|
| Needle Cannula | Stainless steel X5CrNi18-10 |
| Needle hub | Polypropylene in accordance with European Pharmacopoeia Chapter 3.1.3. "Polyolefines". |
| Lubricant | Polydimethylsiloxane complying with Eur. Ph. Chapter 3.1.8 "Silicone oil used as lubricant". |
| Adhesive | Epoxy glue |
| Case | Polypropylene |
| Cap | Polypropylene |
| Sealing label | Self-adhesive paper (tamper-proof closure) |

The needles do not contain components made of natural latex.

Materials of the primary package (foil pouch):

| Component | Material |
|------------|---|
| Foil Pouch | <ul style="list-style-type: none">• Inner layer made of 80-μm cast polypropylene (CPP) used to provide a package seal and maintain seal integrity.• Middle layer made of 18-μm aluminum foil used to provide stiffness and act as a barrier from light, moisture and oxygen.• Outside layer of 12-μm polyester (PET) for reverse-print graphics and to provide cosmetic gloss to the pouches. The foil pouch with three ready seals is covered by polyurethane-based ink, colors: white and blue (panthone 633). |

Materials of the secondary package (carton box):

The cartons are made of a 0.018 solid bleach sulfate (SBS) stock.

Materials of the transport package (carton box):

The carton box is made of cardboard.

TRANSPORTATION, STORAGE, USE CONDITIONS

Transportation is carried out by any means of transport.

RADIESSE Lidocaine injectable implant should be transported, stored and used at a temperature of 15°C - 32°C. The manufacturer does not impose any requirements for atmospheric pressure. Do not use if the expiration date has been exceeded. The expiration date is printed on the product labels.

DISPOSAL

In the territory of the Russian Federation used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance with SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements to the maintenance of the territory of urban and rural settlements, to waterbodies, to drinking water and water supply of citizens, to atmosphere air, to soils, to living areas, to exploitation of manufacturing and communal areas, to organization and realization of sanitary and anti-epidemic (preventive) events" as belonging to hazard class B. Unused syringes and needles must be disposed of as Class A hazardous waste.

REQUIREMENTS FOR ENVIRONMENTAL PROTECTION

With proper disposal, there is no harm to the environment.

WARRANTY

Merz North America, Inc. warrants that reasonable care has been exercised in the design and manufacture of this product.

THIS WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES NOT EXPRESSLY SET FORTH HEREIN, WHETHER EXPRESSED OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ITS PARTICULAR PURPOSE.

Handling and storage of this product, as well as factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Merz North America, Inc.'s control directly affect the product and the results obtained from its use. Merz North America, Inc.'s obligation under this warranty is limited to the replacement of this product and Merz North America, Inc. shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly, arising from the use of this product. Merz North America, Inc. neither assumes, nor authorizes any person to assume for Merz North America, Inc., any other or additional liability or responsibility in connection with this product.




MANUFACTURER:

Merz North America, Inc.
4133 Courtney Street, Suite 10,
Franksville, WI 53126, USA

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE TERRITORY OF THE RUSSIAN FEDERATION:

LLC "Merz Pharma",
125124, Moscow, street 3rd Yamskogo Polya 2, bld. 26, 2nd floor, premises VII, room 29;
Tel. +7 (495) 653 8555; e-mail: info@merz.ru

SYMBOLS

| | | | |
|--|--|--|------------------------------------|
|  | Manufacturer |  | Not Made with Natural Rubber Latex |
|  | Article number |  | Consult instructions for use |
|  | Batch number |  | Do not use if package is damaged |
|  | Use-by date: YYYY-MM |  | Do not re-use |
|  | Caution |  | Do not re-sterilize |
|  | Sterilized using ethylene oxide |  | Temperature limit |
|  | Sterilized using steam or dry heat |  | Needle |
|  | CE mark in accordance with Directive 93/42/EEC relating to medical devices. This mark is followed by the notified body number. |  | Medical device contents |

Copyright © 2024 Merz North America, Inc. All rights reserved. MERZ AESTHETICS and BELOTERO are registered trademarks of Merz Pharma GmbH & Co. KGaA. RADIESSE is a registered trademark of Merz North America, Inc.

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